

ADAMS™ A1c Lite

HA-8380V

Please use only consumables, accessories and after-sales parts provided by ARKRAY.
Please contact your distributor for the correct product codes when purchasing these items.

■ After-sales parts and consumables

● Consumables

Product	Specifications
COLUMN UNIT 80	1 per pack
ELUENT 80A	600 mL, 4 per box
ELUENT 80B	600 mL, 2 per box
ELUENT 80CV	600 mL, 2 per box
HEMOLYSIS WASHING SOLUTION Lite H	2 L, 3 per box
CONTROL DILUTION SET 80	DILUENT 80 (250 mL, 1 per box), RECONSTITUENT 80 (15 mL, 1 per box)
CALIBRATOR Lite	CALIBRATOR Lite Low, CALIBRATOR Lite High and CALIBRATOR Lite Diluent (one for each), Standard value sheet, calibration information barcode labels
Printer paper	58 mm (W) × 25 mm (L), 5 rolls per box
Sample cup	500 µL, 1000 per pack
Sample collection device	100 per pack, for Micro Measurement Kit

● Maintenance parts

Product	Specifications
Power cord	220-240VAC
	125VAC
	The type of power cord (supplied) varies depending on the country. For more information, contact your distributor.
Fuses	T4AE250V~, 2 per pack
Bottle cap A with nozzle	Nozzle for ELUENT 80A
Bottle cap B with nozzle	Nozzle for ELUENT 80B
Bottle cap CV with nozzle	Nozzle for ELUENT 80CV
Bottle cap H with nozzle	Nozzle for Hemolysis Washing Solution Lite H
Cell washing kit	For optical cell washing
Aluminium pack	For maintenance when the instrument is not to be used for extended periods of time, 2 per set
Washing solution for tubes	250 mL
Piercing nozzle	1 per pack
Filter and O-ring for elect.	Mesh filters and O-rings, 5 for each
Silicone tube	2 mm (i.d.) × 4 mm (o.d.), 3 m
Silicone tube	3 mm (i.d.) × 6 mm (o.d.), 3 m
Tefzel tube	1 mm (i.d.) × 2 mm (o.d.), 2 m
Joints 1 × 2	M6 flat seal fitting Ø2, 5 per pack
Joints 2 × 3	M6 flat seal fitting Ø3, 5 per pack
Joint set	M6 flat seal fitting Ø2, through type, 5 per pack
Push screws (circle type)	EASYFIT #9001, 5 per pack
Normal rack	For normal measurements (without adapters)
Calibration rack	For HbA1c calibration, 2 per pack

● Maintenance parts

Product	Specifications
Adapter	Ø13, gray, 10 per pack
Adapter	Ø13, orange, 10 per pack
Adapter	Ø13, blue, 10 per pack

■ Options

Product	Specifications
Liquid waste bottle	With fluid level detection sensor
Hand-held barcode reader	Hand held type
Ethernet board	For connecting the instrument to the host computer
Adapter for Sarstedt tube (Ø11)	Sarstedt tube (Ø11 × L66 mm, 2.7 mL)
Adapter for Sarstedt tube (Ø13)	Sarstedt tube (Ø13 × L65 mm, 2.6 mL)
Micro measurement kit	For micro measurement