

SPOTCHEM™ EZ

SP-4430

Please use only consumables, accessories and after-sales parts provided by ARKRAY.
Please contact your distributor for the correct product codes when purchasing these items.

■ Consumable Accessories

Item	Specifications	Qty
Thermal Printer Paper	For Built-in Paper, 58 mm width	1
Serum Sample Tube	100 pieces, blue cap	1
Serum Sample Tube	500 pieces, blue cap	1
Whole Blood Sample Tube	100 pieces, orange cap	1
Whole Blood Sample Tube	500 pieces, orange cap	1
Pipette Tip	100 pieces, blue	1
Pipette Tip	500 pieces, blue	1
Centrifuge Tube (for SPOTCHEM EZ)	100 pieces	1

■ Parts sold separately

Item	Specifications	Qty
AC Adapter		1
Cleaning Wire	Nozzle cleaning	1
Cleaner Set	Brush, 5 cotton swabs	1
Power cord *	Rating: 125V 13A (A Type Plug)	1
Power cord *	Rating: 250V 10A (C Type Plug)	1
Pipette for solution	For 3 mL	1
Sample Rack	For 30 sample tubes	1
Nozzle Set (EZ)	Replacement Nozzle, O-ring	1
Wrench set for nozzle replacement		1
Tip Waste Case	2 pieces	1
Protective Cover	2 pieces	1

* Choose the power cord that suits your situation.

■ Parts optional

Item	Specifications	Qty
Hand-held barcode reader	Hand held type	1
Portable Centrifuge		1
RS232C Connection Cable	Modular 6P, D-SUB9P, Inch thread	1
Carrying Case	Common Specifications for SPOTCHEM series	1