

Aution|Daten

AE-4070

Please use only consumables, accessories and after-sales parts provided by ARKRAY.
Please contact your distributor for the correct product codes when purchasing these items.

■ Consumables

Product	Specifications
Thermal recording paper	58 mm width, 5 rolls

■ Maintenance parts

Product	Specifications
AC adapter	
Power cord *	125V 7A (A type plug)
Power cord *	250V 2.5A (C type plug)
Test strip tray (AE-4070)	Including suction ports
Waste box (AE-4070)	
Carrying arm (AE-4070)	
Feed lever (AE-4070)	
Test strip feed tray (144)	
Check strip set	White check strips (2) Gray check strips (2)

* Choose the power cord that suits your situation.

■ Options

Product	Specifications
Hand-held barcode reader (AE-4070)	D-SUB9P
Ethernet board (AE-4070)	