

Choose the Model to Fit Your Lifestyle.

GLUCOCARD **W**

GT-7210

Easy to use.

For those who want simplicity

 Blood volume 0.5µL
 Measuring time 7 seconds






GLUCOCARD **S**

GT-7110

Fast and Easy.

For active people who want speed and ease of use

 Blood volume 0.3µL
 Measuring time 5 seconds
 Choice of test strip packaging








GLUCOCARD **Sonyx**

GT-7120

For any situation.

For those who want multiple functions

 Blood volume 0.3µL
 Measuring time 5 seconds
 Choice of test strip packaging
 Test strip guide light & LCD backlight for use in the dark
 Wireless communication function using Bluetooth 4.0.



Features

Equipped with test strip disposal lever

Safe use with no blood contact.



Simple data management using USB communication

Simple blood glucose management.




Flag setting function

Assistance in understanding of measurements.



Accuracy

Compliant with EN ISO 15197:2015



Haematocrit correction function (20-70%)

Peace of mind for patients with polycythemia or anemia.

GLUCOCARD S Test Strips Choice of Packaging Style.



Specifications	GLUCOCARD W	GLUCOCARD S	GLUCOCARD S onyx
Enzyme	GOD	FAD - GDH	
Measurement units	mg/dL or mmol/L (non-switchable)		
Test range	10 - 600 mg/dL, 0.6-33.3mmol/L		
Battery	3 V lithium battery (CR2032) x 1		3V : AAA 1.5V x 2
Memory capacity	500 test results		
Operating environment	Temperature: 8 - 40°C (46 - 104°F) Humidity: 20 - 80 % RH (No condensation)		
Dimensions	L84×W50×H17.6 mm (L3.3×W 2.0×H0.7 inches)		L94×W57×H18.5 mm (L3.7×W 2.2×H0.7 inches)
Weight	Approx. 47 g (1.66 ounces) with battery		Approx. 70 g (2.47 ounces) with battery

*The unit (mg/dL or mmol/L) used with the meter varies in different countries.
You cannot change the preset unit.

In 2017, GLUCOCARD S onyx won the Good Design Award.

Evaluation from the judge

We highly appreciate the efforts to gain trust with the surface texture, design features and UI design and help ordinary people approach unfamiliar medical devices in a more convenient and familiar way.

Yoo's Vision

Its simple yet functional design was aimed to meet ARKRAY's brand identity as well as the best usability since any errors of medical products cannot be tolerated. It comes with the best grip feeling with rounded form and details that offer best usability. To meet premium brand image, its subtle metallic matt black gives understated yet premium look and feel.



Copyright holder : (C)JDP
Site name : GOOD DESIGN AWARD
URL : <http://www.g-mark.org>

YEONGKYU YOO

The founder and creative director of cloudandco, is a highly experienced designer with more than 20 years of successful achievements. While working at Microsoft, he contributed to design work on the HoloLens and Xbox. Before that, he was Creative Director of iriver, and Senior Industrial Designer at NIKE WHQ.

Legal manufacturer

arkray factory, inc.

1480 Koji Konan-cho, Koka, Shiga-ken 520-3306, JAPAN

European representative

arkray europe, b.v.

Prof. J.H. Bavincklaan 2, 1183 AT Amstelveen, THE NETHERLANDS

<https://www.arkray.eu/english/>

*Designs and specifications may be changed without prior notice.

arkray global business, inc.

Yousuien-nai, 59 Gansuin-cho, Kamigyo-ku,

Kyoto 602-0008, JAPAN

TEL +81-75-662-8979 FAX +81-75-431-1202